

Product Name:	ER301WWW/0
Description:	<p>Corepoint Scientific™ General Purpose Hazardous Location (Explosion Proof) Refrigerator 30 Cu. Ft. Capacity</p> <p>One solid swing door, right hinged, not reversible</p> <p>Mechanical thermostat sealed in UL listed explosion proof box</p> <p>Adjustable operating temperature range: 1°C – 10°C</p> <p>Five adjustable shelves</p> <p>Keyed door lock</p> <p>Magnetic door gasket for positive seal</p> <p>Stainless steel door handle</p> <p>Leveling legs</p> <p>White textured cabinet and door</p> <p>High density urethane foam cabinet and door insulation</p> <p>CFC-free refrigerant (R404a), foam insulation and packaging</p> <p>Exterior dimensions: 35 ¾"W x 35"D x 77 ½"H</p> <p>Cycle defrost</p> <p>Not equipped with line cord and plug; must be hardwired in conduit to 120 VAC power source</p> <p>115V, 60Hz, 6 Amps, 1/3 HP</p> <p>UL/C-UL listed</p> <p>One year parts and labor warranty</p> <p>Meet NFPA and OSHA guidelines for hazardous locations, Class 1, Division II, Group C &amp; D.</p> <p>These units are not intended for use in general clinical or research laboratories.</p> <p>Shipping weight: 393 lbs.</p>