





COREPOINT SCIENTIFIC

Interior
FF201WWW/0MHC

Front View



Side View



115VAC/60Hz/1PH
NEMA-5-15P

Cu. Ft	Defrost	Door	Int Door	Shelves	Interior Dim.			Refrigerant	H.P.	Amps	Weight
					W"	D"	H"				
20	Manual	1 Solid	0	6 Fixed	25 ⁷ / ₈	21 ³ / ₄	65	R290	¹ / ₂	6	210 lbs.

Note: Interior dimensions are estimates and subject to change without notice