

|               |  |
|---------------|--|
| Product Name: | LR301WWW/0   |
|               |  |
| Description:  | <p>General Purpose Refrigerator<br/>30 Cu. Ft. Capacity<br/>One solid swing door, right hinged, not reversible<br/>Mechanical thermostat<br/>Temperature range: 0°C to 10°C<br/>Five adjustable shelves<br/>Keyed door lock<br/>Magnetic door gasket for positive seal<br/>Stainless steel door handle<br/>Leveling legs<br/>Probe access port (3/8")<br/>Powder coat paint exterior<br/>High density urethane foam cabinet and door insulation<br/>CFC-free refrigerant (R134a), foam insulation and packaging<br/>Two year parts and labor warranty, plus an additional three year compressor parts warranty<br/>Exterior dimensions: 35 ¾"W x 35" D x 72 ½"H<br/>Auto defrost<br/>115V, 60Hz, 6.8 amps, 1/3 HP<br/>UL/C-UL listed<br/>Shipping weight: 396 lbs.</p> |