





# COREPOINT SCIENTIFIC

## Interior NSWDR492SSS/0

### Front View



### Side View



115VAC/60Hz/1PH  
NEMA-5-15P

| Cu. Ft | Defrost | Door    | Int Door | Shelves | Interior Dim. |        |        | Refrigerant | H.P. | Amps | Weight   |
|--------|---------|---------|----------|---------|---------------|--------|--------|-------------|------|------|----------|
|        |         |         |          |         | W"            | D"     | H"     |             |      |      |          |
| 49     | Auto    | 2 Solid | 0        | 10 Adj  | 49            | 24 5/8 | 52 1/2 | R290        | 3/4  | 9    | 450 lbs. |

Note: Interior dimensions are estimates and subject to change without notice