

Product Name:	PF021WWW/0CAD
Description:	<p>Corepoint Scientific™ White Diamond Series Countertop Controlled Auto Defrost Freezer, Freestanding</p> <p>1 Cu. Ft. Capacity</p> <p>One solid swing door, right hinged</p> <p>Microprocessor temperature controller</p> <p>Digital temperature display with sample simulated probe</p> <p>Audible and visual high and low temperature alarms</p> <p>Remote alarm contacts</p> <p>Operating temperature range: 5°F – -13°F (-15°C – -25°C)</p> <p>Defrost:</p> <p>Auto defrost cycle every 6 hours (at Set Point -24°C)</p> <p>Max air temp in empty freezer: -15°C;</p> <p>Max air temp in loaded freezer: -16°C.</p> <p>Max sample temp in empty freezer: -18°C;</p> <p>Max sample temp in loaded freezer: -20°C.</p> <p>Refrigerator/Freezer freestanding thermometer with 3 year certificate of calibration, with min/max memory, °F/°C switchable</p> <p>2 adjustable shelves</p> <p>Vaccine storage power cord warning label</p> <p>Keyed door lock</p> <p>Magnetic door gasket for positive seal</p> <p>Leveling legs</p> <p>Probe access port (3/8")</p> <p>Powder coat exterior</p> <p>High density urethane foam cabinet and door insulation</p> <p>HFC-free refrigerant (R600a), foam insulation and packaging</p> <p>Pharmacy refrigerator/freezer toolkit and temperature logs</p> <p>Three year parts and labor warranty, plus an additional two year compressor parts only warranty</p> <p>Pyxis®, Omnicell® and AcuDose RX® compatible</p> <p>Exterior dimensions: 23 ¾" W x 24" D x 20 3/8" H</p> <p><i>Note: This unit must have 4" clearance on sides and back for adequate ventilation</i></p> <p>Controlled Auto Defrost, Patented Technology</p> <p>115V, 60 Hz, 1 Amp, 1/8 HP</p> <p>ADA compliant</p> <p>UL/C-UL listed</p> <p>Shipping weight: 106 lbs.</p>