

Product Name:	PF051WWW/0M20
Description:	<p>White Diamond Series Undercounter Laboratory Freezer</p> <p>Freestanding (-20°C)</p> <p>4 Cu. Ft. Capacity</p> <p>One solid swing door, right hinged, not reversible</p> <p>Microprocessor temperature controller</p> <p>Digital temperature display</p> <p>Audible and visual high and low temperature alarms</p> <p>2 probes; 1 in air and 1 in sample bottle</p> <p>Remote alarm contacts</p> <p>Adjustable operating temperature range: -15°C – -25°C</p> <p>2 fixed shelves</p> <p>Bottom footed shelf to allow air circulation</p> <p>Back wall evaporator cover</p> <p>Keyed door lock</p> <p>Magnetic door gasket for positive seal</p> <p>Recessed door handle</p> <p>Leveling legs</p> <p>LED interior light</p> <p>Probe access port (3/8")</p> <p>White cabinet and door</p> <p>High density urethane foam cabinet and door insulation</p> <p>HFC-free refrigerant (R600a), foam insulation and packaging</p> <p>Three year parts and labor warranty, plus an additional two year compressor parts only warranty</p> <p>Exterior dimensions: 23 5/8" W x 24 1/4"D x 32 3/16"H</p> <p><i>Note: This unit must have 4" clearance on sides and back for adequate ventilation</i></p> <p>Manual defrost</p> <p>115V, 60 Hz, 1.4 Amps, 1/8 HP</p> <p>UL/C-UL listed</p> <p>Shipping weight: 128 lbs.</p>