

Product Name:	PRF092SSS/0
Description:	<p>Corepoint Scientific™ White Diamond Series Refrigerator & Freezer Combination</p> <p>9 Cu. Ft. Capacity (4.5 Ref./4.2 Freezer)</p> <p>Two stainless steel exterior swing doors, right hinged* Microprocessor temperature controller, refrigerator and freezer Digital temperature display with sample simulated probe, refrigerator and freezer</p> <p>Adjustable operating temperature range: Refrigerator: 1°C – 10°C; Freezer:-15°C – -25°C</p> <p>Audible and visual high and low temperature alarms, refrigerator and freezer</p> <p>Remote alarm contacts, refrigerator and freezer</p> <p>2 adjustable/1 fixed shelves refrigerator/2 fixed shelves freezer with guard rails on back in Refrigerator to prevent product from being pushed against back wall</p> <p>Back wall evaporator cover, Refrigerator</p> <p>Forced draft circulation Refrigerator</p> <p>Keyed door locks</p> <p>Magnetic door gasket for positive seal</p> <p>Recessed door handles</p> <p>Leveling legs</p> <p>Stacking kit</p> <p>LED interior lights</p> <p>Probe access port (3/8"), 2 ea. (1 for refrigerator/1 for freezer) Stainless steel exterior</p> <p>HFC-free refrigerant (R600a), foam insulation and packaging Three year parts and labor warranty, plus an additional two year compressor parts warranty</p> <p>Exterior dimensions: 23 ¾" W x 24 ½"D x 66 ¾"H</p> <p>Manual defrost freezer / cycle defrost refrigerator</p> <p>115V, 60 Hz, 1.74 Amps Refrigerator, 1.5 Amps Freezer, 1/8 HP Refrigerator and Freezer</p> <p>UL/C-UL listed</p> <p>Shipping weight: 266 lbs.</p> <p> </p> <p>*Please note this unit requires 2 outlets.</p> <p>*This product is available left hinged, product number PRF092SSSLH/0</p>